

ecosash™

Introducing the intelligent new EcoSash™ II, proven sustainable energy savings

At last you can reduce those high energy costs with the intelligent new EcoSash™ II, proven to significantly reduce energy consumption in real laboratory tests.

EcoSash™ II is available on Aircare's Advanced Air™ EcoAir™ fume cupboards, these units come with smart design features that allow the user to have hands free control of the fume cupboard's ultra pressure sensitive sash, improving safety and allowing for more efficient operation.

The intelligent design of the EcoSash™ II allows the unit to know if there is an operator present, this is the key to EcoSash™ II as it allows the unit to offer considerable savings in "make up" air demands (cooling and heating), reducing direct and indirect energy consumption while reducing the dBa noise level within the laboratory.

To find out more about the EcoSash™ II and other energy saving laboratory equipment; or to organise a demonstration of how the EcoSash™ II can benefit your laboratory, visit www.aircare.net.au to view our full range of specialist laboratory equipment.

